



Product Information Sheet

B129 Bacterial Screening Medium 523

Properties

| | |
|--------------------------------|------------------------|
| Form: | Powder |
| Appearance: | Yellow to Cream Powder |
| Application: | Bacterial Screening |
| Solubility: | Water |
| Typical Working Concentration: | 32.15 g/L |
| Storage Temp: | 2 – 6° C |
| Other Notes: | pH = 5.5 – 6.5 |

Formula (mg/L)

| | |
|-------------------------------------|--------|
| Magnesium Sulfate•7H ₂ O | 150 |
| Potassium Phosphate, Monobasic | 2000 |
| Agar | 8000 |
| Casein, Enzymatic Hydrolysate | 8000 |
| Sucrose | 10,000 |
| Yeast Extract | 4000 |

Application Notes

Plant Tissue Culture Tested

References

Revised 3/2007

PhytoTechnology Laboratories, LLC

P.O. Box 12205; Shawnee Mission, KS 66282-2205
Phone: 1-888-749-8682 or 913-341-5343; Fax: 1-888-449-8682 or 913-341-5442
Web Site: www.phytotechlab.com © 2007 PhytoTechnology Laboratories, LLC