




Product Information Sheet

B130



6-Benzylaminopurine Solution (1 mg/mL)

Synonym: BA; N6-Benzyladenine
CAS: 1214-39-7
Formula: C₁₂H₁₁N₅
Molecular Wt: 225.3

Properties

Form: Liquid
Appearance: Colorless, Clear Liquid
Application: Cytokinin
Solubility: Miscible with Water
Typical Working Concentration:
Storage Temp: 2 to 6° C
Storage Temp of Stock Solution:
Other Notes: Plant Tissue Culture Tested

Application Notes

Please Note: While *PhytoTechnology Laboratories™* tests each lot of this product with two or more plant cell/ tissue culture lines, it is the sole responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

References

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 *PhytoTechnology Laboratories, LLC™*