

## PHYTOTECHNOLOGY LABORATORIES' CULTURE CONTAINER SELECTION GUIDE

Product Number	Product Description	Product Composition	Autoclavable	Sterile	Re-useable	Dimensions in mm (round – d x h) (rectangular – l x w x h)	Approx. Surface Area of Medium (cm <sup>2</sup> )	Number of Containers per T077 Tray	Number of Containers per ft <sup>2</sup>	Recommended Closure
P700	PhytoTech Culture Box™, Invertable Container	PS <sup>1</sup>	N	Y	N	100 x 100 x 105	64.	12	9	Included
C185	PC-120™ Culture Bottle, w/ C070 PhytoCap™	PC/PP <sup>5</sup>	Y	N	Y	64 x 64 x 120 42 mm opening	40	30	25	Included (C070)
C168	Culture Vessel, Clear-Con™	PS	N	Y	N	64 x 110	32.2	20	16	Included
C030	Culture Vessel, 8 oz Glass Jar	Glass	Y	N	Y	67 x 77	35.3	20	16	C029
C031	Culture Vessel, 32 oz Glass Jar	Glass	Y	N	Y	67 x 308 (approx)	35.3	20	16	C029
C209	Culture Vessel, 8 oz – Round, Non-graduated	PP <sup>2</sup>	Y	N	Y	115 (rim) x 40 90 (base) x 40	77	6	4	Included
C093	Culture Vessel, 16 oz – Round, Non-graduated	PP	Y	N	Y	115 (rim) x 80 90 (base) x 80	63.6	9	7	Included
C215	Culture Vessel, 16 oz – Round, Non-graduated	PP	Y	N	Y	115 (rim) x 80 90 (base) x 80	63.6	9	7	Included
C221	Culture Vessel, 32 oz – Round, Graduated	PP	Y	N	Y	115 (rim) x 140 90 (base) x 140	63.6	9	7	Included
C532	Culture Vessel, 32 oz – Round, Non-Graduated	PP	Y	N	Y	145 (rim) x 80 122 (base) x 80	116.9	4	4	Included
C032	Culture Bottle, 32 oz, Square Milk Bottle	Glass	Y	N	Y	215 x 84 x 84	144	4	4	S982
E052	ECO 2 Culture Container – Rectangular, XXL Filters	CL - PP <sup>3</sup>	Y	N	Y	125 x 65 x 80	81.3	9	8.5	Included
E048	ECO 2 Culture Container – Rectangular, XL Filters	CL - PP	Y	N	Y	125 x 65 x 80	81.3	9	8.5	Included
E062	ECO 2 Culture Container – Round, XXL Filters	CL - PP	Y	N	Y	100 x 80	78.5	8	6	Included
E058	ECO 2 Culture Container – Round, XL Filters	CL - PP	Y	N	Y	100 x 80	78.5	8	6	Included
C906	Culture Vessel, Magenta GA7	PC/PP <sup>5</sup>	Y	N	Y	77 x 77 x 97	59.3	20	16	Included
C905	Culture Vessel, Magenta GA7-3	PC/PP	Y	N	Y	77 x 77 x 77	59.3	20	16	Included
C930	Culture Tubes, Round Bottom	Glass	Y	N	Y	25 x 150	4.9	—	—	C069, C909, C945
C935	Culture Tubes, Flat Bottom	Glass	Y	N	Y	25 x 95	4.9	—	—	C069, C909, C945
C915	Culture Bag, Star*Pac®, 2 Chamber	HDPE <sup>4</sup>	Y	N	N	6" x 9"	—	—	—	—
C916	Culture Bag, Star*Pac®, 5 Chamber	HDPE	Y	N	N	6" x 9"	—	—	—	—
C918	Culture Bag, Star*Pac®, Single Chamber	HDPE	Y	N	N	6" x 18"	—	—	—	—
C919	Culture Bag, Star*Pac®, Single Chamber	HDPE	Y	N	N	27" x 30"	—	—	—	—

Rev. 2/07

1. PS = Polystyrene
2. PP = Polypropylene
3. CL-PP = Clarified Polypropylene

4. HDPE = High Density Polyethylene
5. PC/PP = Polycarbonate containers w/ Polypropylene lids
6. PS/PP = Polystyrene containers w/ Polypropylene lids