



Product Information Sheet

C195 Coconut Water

Properties

Form:	Liquid
Appearance:	Clear, Colorless to Slight Yellow Tint
Application:	Orchid Media Supplement
Solubility:	Miscible with Water
Typical Working Concentration:	25 to 50 milliliters per liter
Storage Temp:	2 to 6° C
Storage Temp of Stock Solution:	
Other Notes:	Plant Tissue Culture Tested Sterile Filtered

Application Notes

References

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™