



PhytoTechnology Laboratories Introduces the ***PhytoCap***[™]

(Product Number C070)

A New & Improved Alternative To B-Caps



PhytoCap[™] Features and Benefits

- *PhytoCaps*[™] are manufactured with greater temperature compatible plastic which allows for higher temperature tolerance during sterilization.
- The design incorporates durability, stability, and functionality while providing a new level of ergonomic comfort.
- Designed to provide a more secure fit, the *PhytoCap*[™] will fit all standard glass baby food jars.
- The *PhytoCap*[™] is the solution for a cost effective, high quality, baby food jar closure for plant tissue culture.

Orders for the new *PhytoCap*[™] can be made by calling or faxing *PhytoTechnology Laboratories* or your local distributor. Orders can also be placed online at www.phytotechlab.com.



C070	CLOSURE, <i>PHYTOCAP</i>[™] Polypropylene Autoclavable Fits all standard glass baby food jars
------	---

PhytoTechnology Laboratories

PO Box 13481, Shawnee Mission, KS 66282-3481

Phone: 913-341-5343 or 1-888-749-8682

Fax: 913-341-5442 or 1-888-449-8682

www.phytotechlab.com

